

## The Dräger Movita Cable Management System

The integrated Movita Cable Management System will help you organize your workspace and improve your workflow. Developed for use in both OR and ICU environments, the Movita Cable Management System was designed to meet the needs of surgeons, anaesthetists, nurses and cleaning staff alike.



### BETTER WORKFLOW.

By streamlining your work processes, such as cleaning and maintenance, this system will save valuable time and will help you reducing the risk of cable contamination and connection errors.

### CABLE MANAGEMENT FOR IMPROVED ORGANISATION OF WORKPLACES

Ceiling supply units are designed to improve operating room and ICU logistics. Soon after their introduction, the ceiling supply design became the standard worldwide. The Movita Cable Management System now takes this idea a step further by adding integrated cable channels over the entire length of the Movita column. The system stores seldom used cables and hoses in a closed channel, well protected from dust and other potential sources of contamination.


### Streamline your workspace

- Reduce cleaning and disinfection times by storing cables and hoses inside closed cable channels
- Improve workspace organisation by keeping seldom used cables out of the working area
- Reduce the risk of disconnection errors
- Simplify recognition by separating seldom used cables from frequently changed hoses and cables.
- Optimize workflow because of well prepared workplaces
- Reduce the risk of contamination

For maximum cable management capacity, dual cable channels can be installed.


### UPGRADEABILITY

Dräger can upgrade your existing Movita column with the Movita Cable Management System on site with a minimum of downtime. Each Movita column can be equipped with two Movita Cable Management systems of the same length.



Movita with Cable Management System

**The difference is clear**



Movita **without** cable management system



Movita **with** cable management system

**TECHNICAL DATA**

One Movita Cable Management Unit can handle up to four 14 mm gas hoses, as well as ten 7 mm electrical cables. Movita Cable Management is available for the following column lengths:

**Channel length**

500 mm
750 mm
1000 mm
1250 mm
1500 mm

**HEADQUARTERS**

Drägerwerk AG & Co. KGaA  
 Moislinger Allee 53–55  
 23558 Lübeck, Germany

[www.draeger.com](http://www.draeger.com)

**REGION EUROPE CENTRAL AND EUROPE NORTH**

Dräger Medical GmbH  
 Moislinger Allee 53–55  
 23558 Lübeck, Germany  
 Tel +49 451 882 0  
 Fax +49 451 882 2080  
[info@draeger.com](mailto:info@draeger.com)

**REGION MIDDLE EAST, AFRICA, CENTRAL AND SOUTH AMERICA**

Dräger Medical GmbH  
 Branch Office Dubai  
 Dubai Healthcare City, P.O. Box 505108  
 Dubai, United Arab Emirates  
 Tel + 971 436 24 762  
 Fax + 971 436 24 761  
[contactuae@draeger.com](mailto:contactuae@draeger.com)

**REGION NORTH AMERICA**

Draeger Medical, Inc.  
 3135 Quarry Road  
 Telford, PA 18969-1042, USA  
 Tel +1 215 721 5400  
 Toll-free +1 800 437 2437  
 Fax +1 215 723 5935  
[info.usa@draeger.com](mailto:info.usa@draeger.com)

**REGION EUROPE SOUTH**

Dräger Médical S.A.S.  
 Parc de Haute  
 Technologie d'Antony 2  
 25, rue Georges Besse  
 92182 Antony Cedex, France  
 Tel +33 1 46 11 56 00  
 Fax +33 1 40 96 97 20  
[dlimfr-contact@draeger.com](mailto:dlimfr-contact@draeger.com)

**REGION ASIA / PACIFIC**

Draeger Medical  
 South East Asia Pte Ltd  
 25 International Business Park  
 #04-27/29 German Centre  
 Singapore 609916, Singapore  
 Tel +65 6572 4388  
 Fax +65 6572 4399  
[asia.pacific@draeger.com](mailto:asia.pacific@draeger.com)

**Manufacturer:**

Dräger Medical GmbH  
 23542 Lübeck, Germany  
 The quality management system at Dräger Medical GmbH is certified according to ISO 13485, ISO 9001 and Annex II.3 of Directive 93/42/EEC (Medical devices).